

MR202

Electron Microscopy in Materials Characterization

General Introduction, Resolution, Limits on resolution, Lens aberrations

Introduction to SPM/SEM, Electron Optics - Electron Guns and Lenses, Probe diameter and probe current

Electron-Specimen Interactions, Interaction volume, elastic and inelastic scattering





Basics of SEM imaging, Imaging modes, Detectors, Image contrast, Image processing

XEDS and WDS - Principles and practice, Basics

Case studies in Materials Science - Imaging and Analysis

Newer Techniques - EBSD, LVSEM, ESEM

Sample Preparation and a special note about digital imaging/processing



Microscope

Used to 'see' objects not visible to the human eye

Eye can 'resolve' objects $\sim 0.1\text{mm}$ apart

For anything closer, we need a means of magnifying

Note : **BIG** difference between 'seeing' and 'resolving'

Seeing a car approaching (from its headlights)

Resolving the two headlights as separate sources of light

Optical System - Components


Source of Radiation - Visible-light

System of lenses and apertures

Any system used to form an image uses lenses and apertures that have a certain dimension

Limit to resolution arises from the phenomenon of ***diffraction***

Point object is not mapped on to a point - spread out



Diffraction from a single slit

$$\text{Intensity} \sim (\sin(x)/x)^2$$

A big maxima surrounded by smaller maxima

Rayleigh Criterion

Profiles from two adjacent point will overlap

To be able to resolve two points as distinct

$$\theta_R = \sin^{-1}(1.22\lambda/d)$$

This is the *diffraction-limited* resolution limit


To increase resolution

Large d

Small λ

Rayleigh Criterion: a picture

The lens will focus the light to a fuzzy DOT rather than a true point.



Rayleigh Criterion: an example

- Consider the (ideal) resolving ability of the eye
- Estimate D , the diameter of the pupil
- Use $\lambda = 550 \text{ nm}$ (middle of visible spectrum)
- Now calculate the minimum angle the eye can resolve.
- Now calculate how far apart two points of light can be if they are 5 meters away.

Rayleigh Criterion: an example

- with $D = 5 \text{ mm}$ and $\lambda = 550 \text{ nm}$,
$$\theta_{\text{limit}} = \sin^{-1}(1.22 \times 5.5 \times 10^{-7} \text{ m} / .005 \text{ m})$$
$$= 7.7 \times 10^{-3} \text{ degrees}$$
$$= .46 \text{ arc minutes}$$

so $x/L = \tan(\theta_{\text{limit}})$, and

$$x = 5 \text{ m} * \tan(7.7 \times 10^{-3} \text{ degrees}) = .67 \text{ mm}$$


Limits on Resolution

What diameter of telescope would you require to read the numbers on a license plate from a spy satellite?


Limits on Resolution: Telescope

- need to resolve an “x” size of about 1 cm
- “s” is on order of 100 miles or 150 km
- θ_{limit} then must be (in radians)
 $= 1 \text{ cm} / 150 \text{ km} = 7 \times 10^{-8}$
- $\theta_{\text{limit}} = 1.22 \times 5.5 \times 10^{-7} \text{ m} / D = 7 \times 10^{-8}$
 so $D = 10 \text{ m}$ (Hubble has a 2.4 m diameter)

Lens Aberrations - Other limiting factors for resolution



Astigmatism



Electrons as Waves - The particle-wave duality

Electron moving with a velocity ' v ' has a wavelength associated with it

$$\lambda = h/mv \sim 12.247/E(\text{kV})$$


Typical wavelengths


E	λ
100 kV	0.037 Å
200 kV	0.025 Å

Basics of SEM Operation

- Electron gun produces a 'beam'
 - *Thermionic/Field-emission guns*
- Produce a 'tight' spot on the specimen surface
 - *Condenser and Objective lenses*
- Scanning coils 'raster' the beam across the specimen
 - *Size of scan -> Magnification*
- Electron-specimen interactions
 - *Produces a wide variety of signals*
- Detectors to collect the signal
 - *Different detectors for different signals*

Magnification








Integrated circuit bond pads ---2.5kV--- 500X ---SEI Mode




Leaf Cross Section ---10KV--- 100X--- SEI Mode




Some Examples of SEM Images



Electron Gun



Thermionic Emission



$$J_c = A_c T^2 \exp(-E_w/kT)$$

Richardson Equation

$E_w = 4.5 \text{ eV}$ for W, at
 $T = 2700 \text{ K}$


$J = 3.4 \text{ A/cm}^2$

Use 'thermal' energy to excite electrons from a metal

Field Emission


Tip with small radius of curvature

Very high electric fields




Tungsten Hairpin Filament

Work Function = 4.5 eV






7





LaB₆
Filament

Work Function = 2.5 eV

7
7
m
on




Field Emission Gun




Comparison of Electron Sources
(Brightness = Current/area/solid angle)


Source	Brightness (A/cm ² sr)	Source Size	Lifetime (h)	Vacuum Level
Tungsten filament	10 ⁵	30-100 μm	40-100	10 ⁻⁵ Torr
LaB ₆ <100> (2.5 eV)	10 ⁵	5-50 μm	200- 1000	10 ⁻⁷ Torr
Cold FE	10 ⁸	< 5nm	> 1000	10 ⁻¹⁰ Torr
Thermal FE	10 ⁸	< 5 nm	> 1000	10 ⁻⁹ Torr
Schottky	10 ⁸	15-30 nm	> 1000	




Effect of Bias
on the filament
emission



This is what happens when
you turn the filament knob






Comparison of Electron Sources

Note that a very small FE probe carries a much larger current compared to a W filament.

Lenses in the Electron Microscope



The 'Wehnelt' is the first lens in the SEM


Two sets of lenses in the SEM

Larger current through lens \Rightarrow Smaller Focal Length

Condenser lenses - condense the beam to give smaller beam size


Objective lens - final 'probe-forming' lens in the SEM

Larger currents pass through this lens (compared to the condenser lens) - needs to be cooled.



Large working distance

Two types of
Objective Lenses



Very small working
distance

(a) Pin-hole lens - variable working distance, no size limitation of sample, good depth of field

(b) Small focal length => High Resolution

Three important parameters for image formation

Probe Size -

Smaller the probe size, higher the resolution

Smaller the probe size, lower the brightness

Probe Current/Brightness -

Imaging requires 'large' currents

Convergence Angle -

Depth of field increases as angle decreases

Aberrations decrease as angle decreases

Brightness decreases as angle decreases


These three parameters are NOT independent of each other

Probe Size d_p - diameter of beam as it falls on the specimen


Brightness β - Current density/solid angle

Probe convergence α_p - angle made by the cone of electrons with the central axis at the sample surface


Weak
condenser
lens



Strong
condenser
lens



Small
working
distance



Large
working
distance

